


SUPPORTIVE PESSARY SIZE CHART

DIAMETER MEASURED PESSARIES – **ORDER BY SIZE**
AS SHOWN IN THE YELLOW AREA



www.bypmedical.com
Lada sin costo: 01800-624-4175



Ring with Support • KPBS

SIZE	Inch.	mm.
0	1 3/4"	44
1	2"	51
2	2 1/4"	57
3	2 1/2"	64
4	2 3/4"	70
5	3"	76
6	3 1/4"	83
7	3 1/2"	89
8	3 3/4"	95
9	4"	102
10	4 1/4"	108
11	4 1/2"	114
12	4 3/4"	121
13	5"	127

Incontinence Ring • KPON

SIZE	Inch.	mm.
0	1 3/4"	44
1	2"	51
2	2 1/4"	57
3	2 1/2"	64
4	2 3/4"	70
5	3"	76
6	3 1/4"	83
7	3 1/2"	89
8	3 3/4"	95
9	4"	102
10	4 1/4"	108
11	4 1/2"	114
12	4 3/4"	121
13	5"	127

Inflatoball • KPINF

SIZE	Inch.	mm.
S	2"	51
M	2 1/4"	57
L	2 1/2"	64
XL	2 3/4"	70

Donut • KPDO

SIZE in Inches		mm.
2"		51
2 1/4"		57
2 1/2"		64
2 3/4"		70
3"		76
3 1/4"		83
3 1/2"		89
3 3/4"		95

Incontinence Dish • KPCOND


SIZE in mm.	
55	
60	
65	
70	
75	
80	
85	

Flexible Gellhorn • KPGE


SIZE in Inches		mm.
1 1/2"		38
1 3/4"		44
2"		51
2 1/4"		57
2 1/2"		64
2 3/4"		70
3"		76
3 1/4"		83
3 1/2"		89
3 3/4"		95

Shaatz • KPSS

SIZE in Inches		mm.
1 1/2"		38
1 3/4"		44
2"		51
2 1/4"		57
2 1/2"		64
2 3/4"		70
3"		76
3 1/4"		83
3 1/2"		89
3 3/4"		95



* Remember for sizing: The incontinence knob available for this pessary adds 1/2 inch (13 mm) to the diameter.



All pessaries may be ± 1%



CUBE KPEC

SIZE	Inch.	mm.
0	1"	25
1	1 3/16"	30
2	1 3/8"	35
3	1 1/2"	38
4	1 5/8"	41
5	1 3/4"	44
6	2"	51
7	2 1/4"	57



SIZE	In.(A)	In.(B)	mm.(A)	mm.(B)
0	3 1/8"	2"	79	51
1	3 1/4"	2 1/8"	83	54
2	3 1/2"	2 1/4"	89	57
3	3 3/4"	2 3/8"	95	60
4	4 1/4"	2 1/2"	108	64
5	4 1/2"	2 5/8"	114	67
6	4 3/4"	2 3/4"	121	70
7	5"	2 7/8"	127	73
8	5 1/2"	3"	140	76
9	5 3/4"	3 1/8"	146	79



SIZE	In.(A)	In.(B)	In.(C)	mm.(A)	mm.(B)	mm.(C)
0	1 1/16"	1 1/2"	1 3/8"	27	38	35
1	1 1/8"	1 3/4"	1 1/2"	29	44	38
2	1 3/8"	2"	1 5/8"	35	51	41
3	1 1/2"	2 1/8"	1 3/4"	38	54	44
4	1 5/8"	2 1/4"	1 7/8"	41	57	48
5	1 3/4"	2 3/8"	2"	44	60	51
6	1 7/8"	2 1/2"	2 1/8"	48	64	54
7	2"	2 7/8"	2 1/4"	51	73	57
8	2 1/8"	3"	2 3/8"	54	76	60
9	2 1/4"	3 1/8"	2 1/2"	57	79	64
10	2 3/8"	3 1/4"	2 5/8"	60	83	67



SIZE	In.(A)	In.(B)	mm.(A)	mm.(B)
0	2 3/4"	1 3/4"	70	44
1	3"	1 7/8"	76	48
2	3 1/4"	1 15/16"	83	49
3	3 1/2"	2"	89	51
4	3 5/8"	2 1/8"	92	54
5	3 3/4"	2 1/4"	95	57
6	3 7/8"	2 3/8"	98	60
7	4 1/4"	2 1/2"	108	64
8	4 5/8"	2 5/8"	117	67
9	5"	2 3/4"	127	70



SIZE	In.(A)	In.(B)	In.(C)	mm.(A)	mm.(B)	mm.(C)
0	1"	1 15/16"	1 3/8"	25	49	35
1	1 1/16"	2"	1 1/2"	27	51	38
2	1 1/8"	2 1/8"	1 5/8"	29	54	41
3	1 1/4"	2 3/16"	1 3/4"	32	56	44
4	1 3/8"	2 1/4"	1 7/8"	35	57	48
5	1 7/16"	2 3/8"	2"	37	60	51
6	1 1/2"	2 5/8"	2 1/8"	38	67	54
7	1 5/8"	2 7/8"	2 1/4"	41	73	57
8	1 3/4"	3 1/8"	2 3/8"	44	79	60



SHADED AREA SHOWS MINIMUM ASSORTMENT NECESSARY TO FIT 85+% OF PATIENTS

SEE OTHER SIDE